




EasySampler[®]

Stool Collection Kit


Easy and hygienic stool collections





Open the outer container and take out the sample container. Open the sample container.


Important: Leave the absorbent material in the outer container.



Take out the label and with a pen fill it out with:


- A: your name and date of birth
- B: date and time for sample collection.
- C: if stool sample is from more than 1 bowel movement, write how many. Ex. If one: # of BM's= 1 BM

Put the label on the sample container




Put on disposable gloves.

This is to prevent cross contamination




Lift toilet seat and wipe off the toilet rim surface




Unfold EasySampler. Note, there are Five Flaps, orient the one REAR flap towards the back, with the tape side down.

Recommendation: Before fitting bend the paper after the lines




Tape EasySampler onto the back half of the toilet. Starting with the One REAR FLAP, remove the protection paper from the Tape and adhere EasySampler onto the back half of the toilet rim.



Remove the protection paper from the tape, from the TWO MIDDLE FLAPS and adhere them outside the toilet bowl, as far back as possible.

Important: *EasySampler* should not touch the surface of the water




Remove the protection paper from tape of the TWO FRONT FLAPS and adhere them outside the toilet bowl, as close as possible to the middle flaps.

Important: *Make sure that the EasySampler is in the shape of a bowl. If not, re-orient the flaps for the correct positions. Note that tape can be removed from toilet and re-adhered.*




When oriented correctly, lower the Toilet Seat.




Have bowel movement

Important: *Make sure that the stool does not come into contact with toilet water or urine*




Collect stool sample **using** one of the paper spoons attached to the kit **or using** the spoon attached to the container lid.




Using the two paper spoons, push the stool sample into the SAMPLE container. Add **one or two scoops to the container.** -


Important: *The sample container must not be contaminated at the surface*



Remove gloves. Twist them over the paper spoons and dispose. Wash your hands to prevent contamination from sample.




Remove the EasySampler Flaps from the toilet and flush both EasySampler and the stool down toilet.




Screw the lid onto the sample container to seal the sample

Important: *When using watersaving toilets use full flush.
Repeat it until the amount of water are 10-15 liter.
Tape leftovers on the bowl can easily be removed with rubbing alcohol.*



Place the SAMPLE container into the OUTER container, with absorbent material between the two.
Screw the lid on the outer container and place it in the sealing bag. **Seal the bag and place it back in the plastic bag.**



Place the **stool kit and the cool pack** in the **padded envelope**. **Take or send** the stool sample package to the **NIEHS CRU as soon as possible**.